



Enhanced Oral Assessment

# Meet the Light That's Making Headlines

## VELscope Vx!

- Cordless
- Compact
- Affordable

## World's #1 Oral Disease Discovery Device

- Most powerful aid for early detection of oral mucosal abnormalities
- **Viral, fungal, bacterial infections**
- Inflammation from a variety of causes
- **Cancer & pre-cancer**
- **WHO recommended**, allows practices to aspire to an advanced level of patient care
- **FDA approved** and used by over **15,000** practitioners in **23 countries**

## Set a Shining Example!

- Completely non-invasive
- Affordable for any practice
- Compatible with standard digital cameras (iPod 6<sup>th</sup> Generation) for photo documentation
- Entire exam in about two minutes

**The VELscope can even play a role in *saving lives.***

**Life Scanners Medical Solutions Pvt. Ltd.**

Contact No.: +91 9990 444 646 | Email: [hello@velscopeindia.com](mailto:hello@velscopeindia.com)



# The Power of Simplicity


The custom-designed iPod Touch holder magnetically connects to the VELscope, making attachment a simple process. The dedicated application (available from the Apple® App Store™) automatically configures the device for fluorescence photography once the adapter is connected to the VELscope. When the holder is detached from the VELscope, the iPod touch is instantly optimized for white light photography, allowing for flash or continuous illumination image capture.



## Early Visualization - White Light Versus VELscope VX


### Pre-Clinical Discovery

#### Clinical Appearance Under White Light




- Unremarkable under white light


#### VELscope Visualization



- Oddly shaped
- Clearly delineated border
- Biopsy confirmed: Carcinoma-in-Situ (CIS)
- Associated with the HPV Infection



- Is there evidence of tissue abnormality?
- Should a biopsy be ordered?



- Dramatic loss of fluorescence
- Oddly shaped with clearly delineated border
- Biopsy ordered and confirmed Carcinoma-in-Situ (CIS)



LIFE SCANNERS

Life Scanners Medical Solutions Pvt. Ltd.  
Bungalow 65-66/10, Ln 3, Koregaon Park, Pune 411001

+91 9990 444 646